

## SEM Staining Protocol

1. Fix tissue with 2% Glut/ 2% PFA for tissue, 2-2.5% glut for cells
2. Rinse with buffer until all Glut is gone
3. 1% Osmium for 1hr at RT ( this is Hazordous waste). If there is a lot of collagen do Tannic acid
4. Rinse with PBS – this is hazardous waste
5. 0.45um Syringe filter, 1% Glut / 1% tanic acid in PBS for 1hr at RT
6. Rinse well with PBS
7. 1% Osmium at RT for 1hr
8. Rinse well with PBS
9. 1% Glut / 1% tannic acid in PBS for 1hr at RT
10. Rinse well with PBS
11. 1% Osmium at RT for 1hr
12. Rinse well with PBS
13. 70 % EtOH 15 mins x2
14. 95% EtOH 15 mins x2
15. 100% EtOH 20 mins x2 for cells, 30 mins x2 for tissue.
16. Critical point dry
17. Mount
18. Sputter coat
19. Silver paint

Osmium in small benchtop fridge

Tannic acid on chemical shelves

## **Critically Point Dry**

Place sample in an appropriate holder, DO NOT USE permonox, thermonox or polystyrene.

Sample may shrink up to 60 %

Put some EtOH into chamber so samples don't dry out. Make sure no EtOH on rubber.

Finger tight valves, then tighten fully, tighten again

Turn machine on, turn on bone dry CO<sub>2</sub>,

Open cool to -5C, open fill.

Vent to remove EtOH, close vent, fill and re-cool, open fill.

Rest for 5mins

Flush x2. All valves closed.

Turn heat on – go above critical point indicated by red line. When pressure reaches 1500 it self vents. Let cycle x2.

When finished turn heat off, and vent real slow until pressure reaches 250, then open all the way.

Move each screw a little, to release evenly, until they are all off.

Power off, put screws on loose.

## **Mount**

Need Brush , standard carbon adhesive tabs, acetone, stubs, pen ,forceps.

Label stub and clean with acetone.

Stick on tab, use a razor to remove excess.

Place sample on adhesive tab, press around the edge.

Sputter coat

Vortex silver paint, and coat around the edges